

Anti-Beta-Amyloid (1-40) antibody [SRM14787] (STJ11109517)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	IF-P/IHC-P/DB/ELISA
	Host / Source	Rabbit
	Reactivity	Mouse
PRODUCT PROPERTIES		
	Clonality	Monoclonal
	Conjugation	Unconjugated
	Purification	Affinity purification
	Dilution Range	IF-P, 1:200-1:800 IHC-P, 1:500-1:2000 DB, 1:500-1:1000 ELISA, Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
	Isotype	IgG
	Molecular Weight	Protein Weight 91kDa Observed Mw 100kDa
	Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION		
	Gene ID	351
	Gene Symbol	APP
	UniProt ID	A4_HUMAN
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.