

Anti-CD79b antibody [SRM14818] (STJ11109548)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IF/ICC/IF-P/FC/ELISA
Host / Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:1000-1:4000 IF/CC, 1:200-1:2000 IF-P, 1:200-1:2000 FC, 1:500-1:1000 ELISA, Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Isotype	IgG
Molecular Weight	Protein Weight 123kDa Observed Mw 110kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	15985
Gene Symbol	Cd79b
UniProt ID	CD79B_MOUSE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---