

Anti-CD105/Endoglin antibody [SRM14828] (STJ11109558)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/FC
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:500-1:1000 FC, 1:500-1:1000
Isotype	IgG
Molecular Weight	Protein Weight 73kDa Observed Mw 73kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2022
Gene Symbol	ENG
UniProt ID	EGLN_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---