

Anti-IGF2BP1/IMP1 antibody [SRM14839] (STJ11109569)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	WB/IHC-P/IF/ICC/IP/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat
PRODUCT PROPERTIES	
Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:1000-1:6000 IHC-P, 1:200-1:800 IF/ICC, 1:200-1:800 IP, 0.5 Mu g-4 Mu g antibody for 400 Mu g-600 Mu g extracts of whole cells ELISA, Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on you
Isotype	IgG
Molecular Weight	Protein Weight 70kDa Observed Mw 70kDa
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	10642
Gene Symbol	IGF2BP1
UniProt ID	IF2B1_HUMAN
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.