

## Anti-Syk antibody [SRM14840] (STJ11109570)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IP/IF/ICC/IF-P/FC/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:1000-1:2000 IP, 0.5 Mu g-4 Mu g antibody for 400 Mu g-600 Mu g extracts of whole cells IF/ICC, 1:200-1:800 IF-P, 1:200-1:800 FC (intra) , 1:500-1:1000 ELISA, Recommended starting concentration is 1 Mu g/mL. Please optimize t
Isotype	IgG
Molecular Weight	Protein Weight 81kDa Observed Mw 90kDa
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">6850</a>
Gene Symbol	<a href="#">SYK</a>
UniProt ID	<a href="#">KSYK_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---