

Anti-COX5A antibody [SRM14841] (STJ11109571)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	WB/IF/ICC/IF-P/IHC-P/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:20000-1:80000 IF/ICC, 1:200-1:800 IF-P, 1:200-1:800 IHC-P, 1:300-1:3000 ELISA, Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Isotype	IgG
Molecular Weight	Protein Weight 351kDa Observed Mw 350kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9377
Gene Symbol	COX5A
UniProt ID	COX5A_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---