

Anti-Nanog antibody [SRM14858] (STJ11109588)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ChIP/FC
Host / Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:500-1:1000 ChIP, 5 Mu g antibody for 10 Mu g-15 Mu g of Chromatin FC (intra) , 1:500-1:1000
Isotype	IgG
Molecular Weight	Protein Weight 86kDa Observed Mw 80kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	71950
Gene Symbol	Nanog
UniProt ID	NANOG_MOUSE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---