

Anti-DLL4 antibody [SRM14866] (STJ11109596)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IP/FC/ELISA
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:500-1:1000 IP, 0.5 Mu g-4 Mu g antibody for 400 Mu g-600 Mu g extracts of whole cells FC, 1:500-1:1000 ELISA, Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requiremen
Isotype	IgG
Molecular Weight	Protein Weight 30kDa Observed Mw 30kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	54567
Gene Symbol	DLL4
UniProt ID	DLL4_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---