

Anti-AIF1/BA1 antibody [SRM14872] (STJ11109602)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IF-P
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:2500-1:10000 IF-P, 1:200-1:1000
Isotype	IgG
Molecular Weight	Protein Weight 36kDa Observed Mw 34kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	11629
Gene Symbol	Aif1
UniProt ID	AIF1_MOUSE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---