

Anti-Histone H2B antibody [SRM14876] (STJ11109606)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	WB/IHC-P/ELISA/ChIP
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:3000-1:15000 IHC-P, 1:2000-1:4000 ELISA, Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements. ChIP, 3 µg antibody for 5 µg-10 µg of Chromatin
Isotype	IgG
Molecular Weight	Protein Weight 22kDa Observed Mw 15 kDa/23 kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3017/8339/8343/8344/8346/8347/8349
Gene Symbol	H2BC4.H2BC6.H2BC7.H2BC8.H2BC10.H2BC21
UniProt ID	H2B1C_HUMAN H2B2E_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---