

Anti-VDAC3 antibody [SRM14891] (STJ11109621)

GENERAL INFORMATION		
Product Type	Primary antibodies	
Applications	WB/IF/ICC/IF-P/IHC-P/ELISA	
Host / Source	Rabbit	
Reactivity	Human/Mouse/Rat	
PRODUCT PROPERTIES		
Clonality	Monoclonal	
Conjugation	Unconjugated	
Purification	Affinity purification	
Dilution Range	WB, 1:1000-1:4000 IF/ICC, 1:50-1:200 IF-P, 1:50-1:200 IHC-P, 1:50-1:200 ELISA, Recommended starting concentration is 1 ̑g/mL. Please optimize the concentration based on your specific assay requirements.	
Isotype	IgG	
Molecular Weight	Protein Weight 71kDa Observed Mw 80kDa	
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.	
TARGET INFORMATION		
Gene ID	7419	
Gene Symbol	VDAC3	
UniProt ID	VDAC3_HUMAN	
ADDITIONAL INFORMATION		
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.	