

## Anti-Beta-Amyloid antibody [SRM14899] (STJ11109629)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IF-P/IHC-P/DB/ELISA
Host / Source	Rabbit
Reactivity	Mouse

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IF-P, 1:50-1:200 IHC-P, 1:800-1:3200 DB, 1:500-1:1000 ELISA, Recommended starting concentration is 1 $\mu$ g/mL. Please optimize the concentration based on your specific assay requirements.
Isotype	IgG
Molecular Weight	Protein Weight 139kDa Observed Mw 150kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">351</a>
Gene Symbol	<a href="#">APP</a>
UniProt ID	<a href="#">A4_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---