

Anti-HDAC1 antibody [SRM14909] (STJ11109639)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	WB/IP/IF/ICC/IF-P/IHC-P/ChIP/ELISA
	Host / Source	Rabbit
	Reactivity	Human/Mouse
PRODUCT PROPERTIES		
	Clonality	Monoclonal
	Conjugation	Unconjugated
	Purification	Affinity purification
	Dilution Range	WB, 1:5000-1:25000 IP, 0.5 Mu g-4 Mu g antibody for 200 Mu g-400 Mu g extracts of whole cells IF/ICC, 1:2000-1:8000 IF-P, 1:2000-1:8000 IHC-P, 1:5000-1:20000 ChIP, 2-3 Mu g antibody for 10 Mu g-25 Mu g of Chromatin ELISA, Reco
	Isotype	IgG
	Molecular Weight	Protein Weight 72kDa Observed Mw 130kDa
	Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION		
	Gene ID	3065
	Gene Symbol	HDAC1
	UniProt ID	HDAC1_HUMAN
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.