

Anti-GFAP antibody [SRM14922] (STJ11109652)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IF-P
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IF-P, 1:200-1:1000
Isotype	IgG
Molecular Weight	Protein Weight 34kDa Observed Mw 34kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2670
Gene Symbol	GFAP
UniProt ID	GFAP_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---