

Anti-Decorin antibody [SRM14929] (STJ11109659)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IF-P/IHC-P/ELISA
Host / Source	Rabbit
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:1000-1:2000 IF-P, 1:3360-1:13440 IHC-P, 1:800-1:8000 ELISA, Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Isotype	IgG
Molecular Weight	Protein Weight 9kDa Observed Mw Refer to figures
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1634
Gene Symbol	DCN
UniProt ID	PGS2_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---