

**Anti-SGPL1 antibody [SRM14931] (STJ11109661)****GENERAL INFORMATION**

|               |                            |
|---------------|----------------------------|
| Product Type  | Primary antibodies         |
| Applications  | WB/IF/ICC/IF-P/IHC-P/ELISA |
| Host / Source | Rabbit                     |
| Reactivity    | Human/Mouse/Rat            |

**PRODUCT PROPERTIES**

|                     |                                                                                                                                                                                                                            |
|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                                                                                                                                                 |
| Conjugation         | Unconjugated                                                                                                                                                                                                               |
| Purification        | Affinity purification                                                                                                                                                                                                      |
| Dilution Range      | WB, 1:7500-1:30000<br>IF/ICC, 1:200-1:800<br>IF-P, 1:200-1:800<br>IHC-P, 1:100-1:400<br>ELISA, Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements. |
| Isotype             | IgG                                                                                                                                                                                                                        |
| Molecular Weight    | Protein Weight 1011kDa/112kDa<br>Observed Mw 112kDa                                                                                                                                                                        |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                                                                                                                             |

**TARGET INFORMATION**

|             |                             |
|-------------|-----------------------------|
| Gene ID     | <a href="#">8879</a>        |
| Gene Symbol | <a href="#">SGPL1</a>       |
| UniProt ID  | <a href="#">SGPL1_HUMAN</a> |

**ADDITIONAL INFORMATION**

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|