

Anti-METTL3 antibody [SRM14933] (STJ11109663)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	WB/IHC-P/IP/ChIP/ELISA
	Host / Source	Rabbit
	Reactivity	Human/Mouse/Rat
PRODUCT PROPERTIES		
	Clonality	Monoclonal
	Conjugation	Unconjugated
	Purification	Affinity purification
	Dilution Range	WB, 1:5000-1:10000 IHC-P, 1:200-1:2000 IP, 0.5 Mu g-4 Mu g antibody for 200 Mu g-400 Mu g extracts of whole cells ChIP, 5 Mu g antibody for 10 Mu g-15 Mu g of Chromatin ELISA, Recommended starting concentration is 1 Mu g/mL. Please op
	Isotype	IgG
	Molecular Weight	Protein Weight 114kDa Observed Mw 110-125kDa
	Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION		
	Gene ID	56339
	Gene Symbol	METTL3
	UniProt ID	MTA70_HUMAN
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.