

Anti-SMARCAD1 antibody [SRM14936] (STJ11109666)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	WB/IP/IF/ICC/IHC-P/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:25000-1:50000 IP, 0.5 Mu g-4 Mu g antibody for 200 Mu g-400 Mu g extracts of whole cells IF/ICC, 1:200-1:800 IHC-P, 1:200-1:800 ELISA, Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on y
Isotype	IgG
Molecular Weight	Protein Weight 57kDa Observed Mw 66kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	56916
Gene Symbol	SMARCAD1
UniProt ID	SMRCD_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---