

Anti-Pro-Collagen I/COL1A1 antibody [SRM14940] (STJ11109670)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	WB/IP/IF-P/IHC-P/ELISA
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:13000-1:78000 IP, 0.5 Mu g-4 Mu g antibody for 200 Mu g-400 Mu g extracts of whole cells IF-P, 1:200-1:800 IHC-P, 1:1000-1:10000 ELISA, Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on
Isotype	IgG
Molecular Weight	Protein Weight 14kDa Observed Mw 14kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1277
Gene Symbol	COL1A1
UniProt ID	CO1A1_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---