

Anti-TMEM119 antibody [SRM14947] (STJ11109677)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IF-P/IHC-P/mIHC/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:2000-1:24000 IF-P, 1:500-1:2000 IHC-P, 1:2000-1:8000 mIHC, 1:2000-1:8000 ELISA, Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Isotype	IgG
Molecular Weight	Protein Weight 31kDa Observed Mw 31kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	338773
Gene Symbol	TMEM119
UniProt ID	TM119_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---