

Anti-IRF3 antibody [SRM14968] (STJ11109698)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-P/IF/ICC/IP/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:11000-1:44000 IHC-P, 1:600-1:3000 IF/ICC, 1:100-1:400 IP, 0.5 μ g-4 μ g antibody for 200 μ g-400 μ g extracts of whole cells ELISA, Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on
Isotype	IgG
Molecular Weight	Protein Weight 75 kDa Observed Mw 75 kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3661
Gene Symbol	IRF3
UniProt ID	IRF3_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	-------------------------------------------------------------------------------------------------------------------------