

Anti-ECHS1 antibody [SRM14976] (STJ11109706)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	WB/IHC-P/IF/ICC/ELISA
	Host / Source	Rabbit
	Reactivity	Human/Mouse/Rat
PRODUCT PROPERTIES		
	Clonality	Monoclonal
	Conjugation	Unconjugated
	Purification	Affinity purification
	Dilution Range	WB, 1:2000-1:12000 IHC-P, 1:2000-1:8000 IF/ICC, 1:200-1:800 ELISA, Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
	Isotype	IgG
	Molecular Weight	Protein Weight 44kDa Observed Mw 53kDa
	Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION		
	Gene ID	1892
	Gene Symbol	ECHS1
	UniProt ID	ECHM_HUMAN
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.