

Anti-AIF1/BA1 antibody [SMM14980] (STJ11109710)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	WB/IF-F/IF-P/IHC-P/ELISA
Host / Source	Mouse
Reactivity	Human/Mouse/Rat
PRODUCT PROPERTIES	
Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:1000-1:3000 IF-F, 1:100-1:400 IF-P, 1:100-1:400 IHC-P, 1:10000-1:50000 ELISA, Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Isotype	IgG
Molecular Weight	Protein Weight 45kDa Observed Mw 50kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	11629
Gene Symbol	Aif1
UniProt ID	AIF1_MOUSE
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.