

Anti-Lactate Dehydrogenase B/LDH-B antibody [SRM14993] (STJ11109723)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IF/ICC/IF-P/IHC-P/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:12000-1:48000 IF/CC, 1:200-1:1600 IF-P, 1:200-1:1600 IHC-P, 1:5000-1:20000 ELISA, Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Isotype	IgG
Molecular Weight	Protein Weight 67kDa Observed Mw 67kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3945
Gene Symbol	LDHB
UniProt ID	LDHB_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---