

Anti-VGLUT2 antibody [SRM14994] (STJ11109724)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IF/ICC/IF-F/IF-P/IHC-P/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:3000-1:12000 IF/ICC, 1:200-1:400 IF-F, 1:200-1:400 IF-P, 1:200-1:400 IHC-P, 1:1500-1:6000 ELISA, Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Isotype	IgG
Molecular Weight	Protein Weight 40kDa Observed Mw 40kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	57084
Gene Symbol	SLC17A6
UniProt ID	VGLU2_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---