

Anti-Myelin Basic Protein antibody [SMM14996] (STJ11109726)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	WB/IF-P/IHC-P/ELISA
	Host / Source	Mouse
	Reactivity	Human/Mouse/Rat
PRODUCT PROPERTIES		
	Clonality	Monoclonal
	Conjugation	Unconjugated
	Purification	Affinity purification
	Dilution Range	WB, 1:10000-1:70000 IF-P, 1:200-1:2000 IHC-P, 1:2000-1:10000 ELISA, Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
	Isotype	IgG
	Molecular Weight	Protein Weight 61kDa Observed Mw 55kDa
	Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION		
	Gene ID	4155
	Gene Symbol	MBP
	UniProt ID	MBP_HUMAN
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.