

Anti-Rbm7 antibody [SRM14999] (STJ11109729)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IP/IF-P/IHC-P/ELISA
Host / Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:2500-1:10000 IP, 0.5 Mu g-4 Mu g antibody for 200 Mu g-400 Mu g extracts of whole cells IF-P, 1:200-1:1000 IHC-P, 1:2000-1:8000 ELISA, Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on y
Isotype	IgG
Molecular Weight	Protein Weight 40kDa Observed Mw 49kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	67010
Gene Symbol	Rbm7
UniProt ID	RBM7_MOUSE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---