

Anti-Synaptophysin antibody [SRM15001] (STJ11109731)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IP/IF-F/IF-P/IHC-P/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:5000-1:40000 IP, 0.5 Mu g-4 Mu g antibody for 500 Mu g-700 Mu g extracts of whole cells IF-F, 1:200-1:800 IF-P, 1:200-1:800 IHC-P, 1:500-1:2000 ELISA, Recommended starting concentration is 1 Mu g/mL. Please optimize the conc
Isotype	IgG
Molecular Weight	Protein Weight 120kDa Observed Mw 145kDa
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6855
Gene Symbol	SYP
UniProt ID	SYPH_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---