

Anti-PML antibody [SRM15011] (STJ11109741)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IP/IF/ICC/IHC-P/ELISA
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:2500-1:10000 IP, 0.5 Mu g-4 Mu g antibody for 200 Mu g-400 Mu g extracts of whole cells IF/ICC, 1:200-1:800 IHC-P, 1:5000-1:20000 ELISA, Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on
Isotype	IgG
Molecular Weight	Protein Weight 42kDa Observed Mw 75-100kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5371
Gene Symbol	PML
UniProt ID	PML_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---