

Anti-AMPK alpha 1/2 antibody [SRM15017] (STJ11109747)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	WB/IF/ICC/IHC-P/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:12000-1:50000 IF/CC, 1:100-1:200 IHC-P, 1:5000-1:20000 ELISA, Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Isotype	IgG
Molecular Weight	Protein Weight 36kDa Observed Mw 36kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5562 5563
Gene Symbol	PRKAA1 PRKAA2
UniProt ID	AAPK1_HUMAN AAPK2_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---