

Anti-TIM23 antibody [SRM15022] (STJ11109752)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	WB/IP/IF/ICC/IHC-P/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:5000-1:30000 IP, 0.5 Mu g-4 Mu g antibody for 200 Mu g-400 Mu g extracts of whole cells IF/ICC, 1:200-1:400 IHC-P, 1:500-1:2000 ELISA, Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on y
Isotype	IgG
Molecular Weight	Protein Weight 25kDa Observed Mw 25kDa
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	100287932
Gene Symbol	TIMM23
UniProt ID	TIM23_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---