

Anti-SLC16A3/MCT4 antibody [SRM15023] (STJ11109753)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	WB/IF/ICC/IHC-P/ELISA
	Host / Source	Rabbit
	Reactivity	Human
PRODUCT PROPERTIES		
	Clonality	Monoclonal
	Conjugation	Unconjugated
	Purification	Affinity purification
	Dilution Range	WB, 1:5000-1:10000 IF/ICC, 1:100-1:200 IHC-P, 1:200-1:800 ELISA, Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
	Isotype	IgG
	Molecular Weight	Protein Weight 25kDa Observed Mw 23kDa
	Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION		
	Gene ID	9123
	Gene Symbol	SLC16A3
	UniProt ID	MOT4_HUMAN
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.