

Anti-Fibrillarin antibody [SRM15035] (STJ11109765)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	WB/IP/IF/ICC/IHC-P/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat
PRODUCT PROPERTIES	
Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:25000-1:100000 IP, 0.5 Mu g-4 Mu g antibody for 200 Mu g-400 Mu g extracts of whole cells IF/ICC, 1:200-1:400 IHC-P, 1:1000-1:4000 ELISA, Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based o
Isotype	IgG
Molecular Weight	Protein Weight 74kDa Observed Mw 80kDa
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	2091
Gene Symbol	FBL
UniProt ID	FBRL_HUMAN
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.