

Anti-CDKN1A/p21CIP1 antibody [SRM15042] (STJ11109772)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IP/IF/ICC/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:1000-1:5000 IP, 0.5 Mu g-4 Mu g antibody for 500 Mu g-700 Mu g extracts of whole cells IF/ICC, 1:100-1:400 ELISA, Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requir
Isotype	IgG
Molecular Weight	Protein Weight 85kDa Observed Mw 101kDa
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1026
Gene Symbol	CDKN1A
UniProt ID	CDN1A_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---