

**Anti-Pit1/POU1F1 antibody [SRM15044] (STJ11109774)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	WB/IP/IF-P/IHC-P/ELISA
Host / Source	Rabbit
Reactivity	Mouse/Rat

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:4000-1:16000 IP, 0.5 Mu g-4 Mu g antibody for 200 Mu g-400 Mu g extracts of whole cells IF-P, 1:200-1:800 IHC-P, 1:1000-1:4000 ELISA, Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on yo
Isotype	IgG
Molecular Weight	Protein Weight 435kDa Observed Mw 40-100kDa
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

**TARGET INFORMATION**

Gene ID	<a href="#">5449</a>
Gene Symbol	<a href="#">POU1F1</a>
UniProt ID	<a href="#">PIT1_HUMAN</a>

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---