

Anti-ATF6 antibody [SRM15045] (STJ11109775)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IP/IHC-P/ELISA
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:5000-1:10000 IP, 0.5 Mu g-4 Mu g antibody for 200 Mu g-400 Mu g extracts of whole cells IHC-P, 1:200-1:800 ELISA, Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requir
Isotype	IgG
Molecular Weight	Protein Weight 68kDa Observed Mw 70kDa
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	22926
Gene Symbol	ATF6
UniProt ID	ATF6A_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---