

Anti-STAT3 antibody [SRM15052] (STJ11109782)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IF/ICC/IHC-P/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat/Monkey

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:5000-1:20000 IF/CC, 1:200-1:400 IHC-P, 1:1000-1:4000 ELISA, Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Isotype	IgG
Molecular Weight	Protein Weight 28kDa Observed Mw 17kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6774
Gene Symbol	STAT3
UniProt ID	STAT3_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---