

Anti-IL-6 antibody [SRM15055] (STJ11109785)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IP/IF/ICC/ELISA
Host / Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:1000-1:5000 IP, 0.5 Mu g-4 Mu g antibody for 900 Mu g-1100 Mu g extracts of whole cells IF/ICC, 1:200-1:800 ELISA, Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requi
Isotype	IgG
Molecular Weight	Protein Weight 45kDa Observed Mw 45-50kDa
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	24498
Gene Symbol	IL6
UniProt ID	IL6_RAT

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---