

Anti-CD274 antibody [SRM15057] (STJ11109787)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IF/ICC/IHC-P/FC/ELISA
Host / Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:1000-1:4000 IF/CC, 1:200-1:400 IHC-P, 1:200-1:800 FC, 1:100-1:500 ELISA, Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements. For high-ratio antibody dilu
Isotype	IgG
Molecular Weight	Protein Weight 861kDa Observed Mw 55kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	60533
Gene Symbol	Cd274
UniProt ID	PD1L1_MOUSE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---