

**Anti-hCG beta antibody [SRM15063] (STJ11109793)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	WB/IF-P/IHC-P/ELISA
Host / Source	Rabbit
Reactivity	Human

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:2500-1:5000 IF-P, 1:200-1:400 IHC-P, 1:200-1:800 ELISA, Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Isotype	IgG
Molecular Weight	Protein Weight 109kDa Observed Mw 125kDa/130kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

**TARGET INFORMATION**

Gene ID	<a href="#">1082/93659/94115</a>
Gene Symbol	<a href="#">CGB3.CGB5.CGB8</a>
UniProt ID	<a href="#">CGB3_HUMAN</a>

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	-------------------------------------------------------------------------------------------------------------------------