

Anti-DDDDK-Tag antibody [SRM15108] (STJ11109838)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IF/ICC/IP/FC/ELISA/ChIP/ChIP-seq
Host / Source	Rabbit
Reactivity	Species Independent

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:5000-1:20000 IF/ICC, 1:50-1:200 IP, 0.5 Mu g-4 Mu g antibody for 200 Mu g-400 Mu g extracts of whole cells FC, 1:50-1:200 ELISA, Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your sp
Isotype	IgG
Molecular Weight	Protein Weight 12kDa Observed Mw 12kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---