

Anti-GFP-Tag antibody [SMM15111] {Affinity Gel} (STJ11109841)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	IP
Host / Source	Mouse
Reactivity	Species Independent

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Affinity Gel
Purification	Affinity purification
Dilution Range	IP, 20ul antibody (bead slurry) for 200 Mu g-400 Mu g extracts of whole cells
Isotype	IgG2bk
Molecular Weight	Protein Weight 141kDa Observed Mw 30-35kDa
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---