

Anti-GFP-Tag antibody [SMM15114] (STJ11109844)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IP
Host / Source	Mouse
Reactivity	Species Independent

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:25000-1:100000 IP, 0.5 Mu g-4 Mu g antibody for 100 Mu g-300 Mu g extracts of whole cells
Isotype	IgG2a
Molecular Weight	Protein Weight 141kDa Observed Mw 10kDa
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---