

Anti-Phospho-Tau-T181-T181 antibody [SRM15126] (STJ11109856)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IP/IF-P/IHC-P/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:1000-1:5000 IP, 0.5 Mu g-4 Mu g antibody for 200 Mu g-600 Mu g extracts of whole cells IF-P, 1:50-1:200 IHC-P, 1:50-1:200 ELISA, Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your sp
Isotype	IgG
Molecular Weight	Protein Weight 14kDa Observed Mw 15kDa
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4137
Gene Symbol	MAPT
UniProt ID	TAU_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---