

Anti-Phospho-LATS1-T1079+LATS2-T1041-T1079+T1041 antibody [SRM15143] (STJ11109873)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	WB/IF/ICC/IHC-P/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:1000-1:2000 IF/ICC, 1:200-1:800 IHC-P, 1:50-1:200 ELISA, Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Isotype	IgG
Molecular Weight	Protein Weight 58kDa Observed Mw 58kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9113 26524
Gene Symbol	LATS1 LATS2
UniProt ID	LATS1_HUMAN LATS2_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---