

## Anti-CIRH1A antibody (497-686) (STJ111325)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse

### PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 77kDa Observed Mw: 77kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">84916</a>
Gene Symbol	<a href="#">UTP4</a>
UniProt ID	<a href="#">UTP4_HUMAN</a>
Immunogen Region	497-686
Immunogen Sequence	NWLAASGTSAGVHVYVNVKQL KLHCTVPAYNFPVTAMAIAIP NTNNLVIAHSDQQVFEYSIP DKQYTDWSRTVQKQGFHHLW LQRDTPITHISFHPKRPMMHI LLHDAYMFCIIDKSLPLPND KTLLYNPFPTNESDVIRRR TAHAFKISKIYKPLLFMDLL DERTLVAVERPLDDIIAQLP PPIKPKKFGT
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 497-686 of human CIRH1A (NP_116219.2).

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of various lysates using CIRH1A Rabbit pAb (STJ111325) at 1:1000 dilution.

Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution.

Lysates/proteins: 25 µg per lane.  
Blocking buffer: 3% nonfat dry milk in TBST.  
Detection: ECL Basic Kit  
Exposure time: 5s.