

## Anti-BBS10 antibody (474-723) (STJ111434)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse

### PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:200-1:3000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 81kDa Observed Mw: 80kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">79738</a>
Gene Symbol	<a href="#">BBS10</a>
UniProt ID	<a href="#">BBS10_HUMAN</a>
Immunogen Region	474-723
Immunogen Sequence	AENKDALEKTQTYLKVHNSL VIPDVELETYIPYSTPTLTP TDTFQTVETLTCLSLERNRL TDYYEPLLKNNSTAYSTRGN RIEISYENLQVTNITRKGSM LPVSCKLPNMGTSQSYLSSS MPAGCVLPVGGNFEILLHY LLNYAKKCHQSEETMVMII ANALLGIPKVLKSKTGKYS FPHTYIRAVHALQTNQPLVS SQTGLESVMGKYQLLTSVLQ CLTKILTIDMVTVKRHPQ
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 474-723 of human BBS10 (NP_078961.3).

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of various lysates using BBS10 Rabbit pAb (STJ111434) at 1:1000 dilution.  
Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution.  
Lysates/proteins: 25 µg per lane.  
Blocking buffer: 3% nonfat dry milk in TBST.  
Detection: ECL Basic Kit  
Exposure time: 3s.