

Anti-VBP1 antibody (1-197) (STJ111509)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:1000-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 23kDa Observed Mw: 23kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7411
Gene Symbol	VBP1
UniProt ID	PFD3_HUMAN
Immunogen Region	1-197
Immunogen Sequence	MAAVKDSCGKGEMATGNGRR LHLGIPEAVFVEDVDSFMKQ PGNETADTVLKKLDEQYQKY KFMEINLAQKKRRLKGQIPE IKQTLKILKYMOKKKKSTNS METRFLADNLYCKASVPPT DKVCLWLGANVMLEYDIDEA QALLEKNLSTATKNLDSLEE DLDFLRDQFTTTEVNMARVY NWDVKRRNKDDSTKNKA
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-197 of human VBP1 (NP_003363.1).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Western blot analysis of various lysates using VBP1 Rabbit pAb (STJ111509) at 1:1000 dilution.
 Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution.
 Lysates/proteins: 25 µg per lane.
 Blocking buffer: 3% nonfat dry milk in TBST.
 Detection: ECL Basic Kit
 Exposure time: 1s.