

Anti-MAD2L2 antibody (1-211) (STJ111903)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 24kDa Observed Mw: 25kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10459
Gene Symbol	MAD2L2
UniProt ID	M2L2_HUMAN
Immunogen Region	1-211
Immunogen Sequence	MTTLTRQDLNFGQVVADVLC EFLEVAVHLILYVREVYPVG IFQKRKKYNVPVQMSCHPEL NQYIQDTLHCVKPLLEKNDV EKVVVVILDKEHRPVEKFVF EITQPPLLSISSDSLSSHVE QLLRAFILKISVCDVLDHN PPGCTFTVLVHTREAATRNM EKIQVIKDFPWILADEQDVH MHDPRLIPLKMTSDILKMQ LYVEERAHKGS
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-211 of human MAD2B/MAD2B/MAD2L2 (NP_006332.3).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Western blot analysis of various lysates using MAD2B/MAD2L2 Rabbit pAb (STJ111903) at 1:1000 dilution.
 Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution.
 Lysates/proteins: 25 µg per lane.
 Blocking buffer: 3% nonfat dry milk in TBST.
 Detection: ECL Basic Kit
 Exposure time: 5s.