

## Anti-POMGNT1 antibody (411-660) (STJ111921)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

### PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 $\mu$ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 75kDa Observed Mw: 90kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">55624</a>
Gene Symbol	<a href="#">POMGNT1</a>
UniProt ID	<a href="#">PMGT1_HUMAN</a>
Immunogen Region	411-660
Immunogen Sequence	HLLEEDDSLYCISAWNDQGY EHTAEDPALLYRVETMPGLG WVLRRLSYKEELEPKWPTPE KLWDWDMWMRMPEQRRGREG IIPDVSRSYHFGIVGLNMNG YFHEAYFKKHKFNTVPGVQL RNVDSLKKEAYEVEVHRLLS EAEVLDHSKNPCEDSFLPDT EGHTYVAFIRMEKDDDDFTTW TQLAKCLHIWDLVVRGNHRG LWRLFRKKNHFLVVGVPASP YSVKKPPSVTPIFLEPPPK
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 411-660 of human POMGNT1 (NP_060209.3).

### ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**



Western blot analysis of various lysates using POMGNT1 Rabbit pAb (STJ111921) at 1:1000 dilution.

Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution.

Lysates/proteins: 25  $\mu$ g per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit

Exposure time: 60s.